







Product Information

Product	<u>Men & Women Sizes</u>	<u>Other Information</u>
<p data-bbox="321 384 412 415">Titan</p> 	<p data-bbox="680 489 771 520">Sizes:</p> <ul data-bbox="680 531 792 562" style="list-style-type: none"> • 5 - 17 <p data-bbox="680 615 797 646">Widths:</p> <ul data-bbox="680 657 824 783" style="list-style-type: none"> • Medium • Wide • X - Wide <p data-bbox="659 867 964 898">Colors: Black or Beige</p> <p data-bbox="659 940 873 972">Closure: Velcro</p> <p data-bbox="659 1003 967 1035">Cover: Lycra / Zennon</p>	<p data-bbox="1101 489 1520 678">Cover: Entire shoe covered in Zennon – a breathable spandex like material. Zennon, gently wraps the foot, giving a feeling a wearing bedroom slippers.</p> <p data-bbox="1101 709 1520 961">Specifications: Indoor shoe only. Ability to wash shoes, bag included. Only foot measurements needed to determine size. Recommended for people with corns, calluses, claw feet, hammer or cross toes.</p>
<p data-bbox="298 1098 435 1129">Zennon</p> 	<p data-bbox="680 1203 771 1234">Sizes:</p> <ul data-bbox="680 1245 792 1276" style="list-style-type: none"> • 4 - 15 <p data-bbox="680 1329 797 1360">Widths:</p> <ul data-bbox="680 1371 844 1549" style="list-style-type: none"> • Medium • Wide • X - Wide • XX - Wide <p data-bbox="659 1633 997 1696">Colors: Black, Brown, or Taupe</p> <p data-bbox="659 1738 980 1770">Closure: Velcro or Lace</p> <p data-bbox="659 1801 967 1833">Cover: Lycra / Zennon</p>	<p data-bbox="1101 1203 1520 1623">Specifications: Same shoe pattern as Leather, Suede, Boot and New Yorker. Recommended for people with corns, calluses, and foot deformities like claw feet, hammer, and/or cross toes. Non-specific ball area on sole. Thick toe box for added protection. Loose tongue that will stretch with foot swelling around ankle and/or front of foot</p>


Zennon		
	<p>Sizes:</p> <ul style="list-style-type: none"> • 4 - 15 <p>Widths:</p> <ul style="list-style-type: none"> • Medium • Wide • X - Wide • XX - Wide <p>Colors: Black, Brown, or Taupe</p> <p>Closure: Velcro or Lace</p> <p>Cover: Lycra / Zennon</p>	<p>Specifications: Same shoe pattern as Leather, Suede, Boot and New Yorker. Recommended for people with corns, calluses, and foot deformities like claw feet, hammer, and/or cross toes. Non-specific ball area on sole. Thick toe box for added protection. Loose tongue that will stretch with foot swelling around ankle and/or front of foot.</p>
SPORT I		
	<p>Sizes:</p> <ul style="list-style-type: none"> • 5 - 15 <p>Widths:</p> <ul style="list-style-type: none"> • Medium • Wide • X - Wide <p>Colors: Black or White</p> <p>Closure: Velcro or Lace</p> <p>Cover: Leather</p>	<p>Specifications: Same Shoe pattern as Sport II. Recommended walking shoe. Non-slippery sole. High back heel padding for narrow heel patients. Padded tongue for high arch and high volume feet.</p>



SPORT I		
	<p>Sizes:</p> <ul style="list-style-type: none"> • 5 - 15 <p>Widths:</p> <ul style="list-style-type: none"> • Medium • Wide • X - Wide <p>Colors: Black or White</p> <p>Closure: Velcro or Lace</p> <p>Cover: Leather</p>	<p>Specifications: Same Shoe pattern as Sport II. Recommended walking shoe. Non-slippery sole. High back heel padding for narrow heel patients. Padded tongue for high arch and high volume feet.</p>
SPORT II		
	<p>Sizes:</p> <ul style="list-style-type: none"> • 5 - 15 <p>Widths:</p> <ul style="list-style-type: none"> • Medium • Wide • X - Wide <p>Colors: Black or White</p> <p>Closure: Velcro or Lace</p> <p>Cover: Napa Leather</p>	<p>Specifications: Same Shoe pattern as Sport I. Recommended walking shoe. Non-slippery sole. Padded tongue for high arch and high volume feet.</p>

Leather		
	<p>Sizes:</p> <ul style="list-style-type: none"> • 4 - 15 <p>Widths:</p> <ul style="list-style-type: none"> • Medium • Wide • X - Wide • XX - Wide <p>Colors: Black, White, or Beige</p> <p>Closure: Velcro or Lace</p> <p>Cover: Leather</p>	<p>Specifications: Same shoe pattern as Zennon, Suede, Boot and New Yorker. Recommended for feet with lots of volume. Non-specific ball area on sole. Thick toe box for added protection. Leather provides additional support around foot.</p>
Leather		
	<p>Sizes:</p> <ul style="list-style-type: none"> • 4 - 15 <p>Widths:</p> <ul style="list-style-type: none"> • Medium • Wide • X - Wide • XX - Wide <p>Colors: Black, White, or Beige</p> <p>Closure: Velcro or Lace</p> <p>Cover: Leather</p>	<p>Specifications: Same shoe pattern as Zennon, Suede, Boot and New Yorker. Recommended for feet with lots of volume. Non-specific ball area on sole. Thick toe box for added protection. Leather provides additional support around foot.</p>

Suede		
	<p>Sizes:</p> <ul style="list-style-type: none"> • 4 - 15 <p>Widths:</p> <ul style="list-style-type: none"> • Wide <p>Colors: Dark Brown</p> <p>Closure: Lace</p> <p>Cover: Suede</p>	<p>Specifications: Same shoe pattern as Zennon, Leather, Boot and New Yorker. Recommended for feet with lots of volume. Non-specific ball area on sole. Thick toe box for added protection.</p>
Boot		
	<p>Sizes:</p> <ul style="list-style-type: none"> • 7 - 14 <p>Widths:</p> <ul style="list-style-type: none"> • Wide • X - Wide <p>Colors: Black</p> <p>Closure: Lace</p> <p>Cover: Leather</p>	<p>Specifications: Same shoe pattern as Zennon, Leather, Suede, and New Yorker. Thick toe box for added protection. Hi-top leather provides additional support around foot and ankle.</p>

Jenny		
	<p>Sizes:</p> <ul style="list-style-type: none"> • 5 - 11 <p>Widths:</p> <ul style="list-style-type: none"> • Narrow • Medium • Wide • X - Wide <p>Colors: Black or Taupe</p> <p>Closure: Velcro</p> <p>Cover: Napa Leather or Lycra / Zennon</p>	<p>Specifications: Narrow width ladies dress shoe. Velcro closure allows for some top swelling.</p>
Dawn		
	<p>Sizes:</p> <ul style="list-style-type: none"> • 5 - 11 <p>Widths:</p> <ul style="list-style-type: none"> • Narrow • Medium • Wide <p>Colors: Black</p> <p>Closure: Velcro</p> <p>Cover: Napa Leather</p>	<p>Specifications: Narrow shoe pattern same as Kosmo. Very stable sole structure.</p>

<p style="text-align: center;">Raven</p>		
	<p>Sizes:</p> <ul style="list-style-type: none"> • 5 - 12 <p>Widths:</p> <ul style="list-style-type: none"> • Medium • Wide • X - Wide <p>Colors: Black or Beige</p> <p>Closure: Velcro</p> <p>Cover: Leather & Lycra / Zennon</p>	<p>Specifications: Narrow width shoe. Toe angle creates shallow toe box. Requires larger shoe size to accommodate for narrow shoe style.</p>
<p style="text-align: center;">Kosmo</p>		
	<p>Sizes:</p> <ul style="list-style-type: none"> • 7 - 14 <p>Widths:</p> <ul style="list-style-type: none"> • Wide • X - Wide • XX - Wide <p>Colors: Black</p> <p>Closure: Lace</p> <p>Cover: Napa Leather</p>	<p>Specifications: Narrow shoe pattern same as Dawn. Very stable sole structure.</p>

<p align="center">Gator</p>		
	<p>Sizes:</p> <ul style="list-style-type: none"> • 5 - 13 <p>Widths:</p> <ul style="list-style-type: none"> • Medium • Wide • X - Wide <p>Colors: Black or Brown</p> <p>Closure: Velcro or Lace</p> <p>Cover: Fine Leather</p>	<p>Specifications: Narrow width dress shoe. Top stitching creates shallow toe box.</p>
<p align="center">Gator</p>		
	<p>Sizes:</p> <ul style="list-style-type: none"> • 5 - 13 <p>Widths:</p> <ul style="list-style-type: none"> • Medium • Wide • X - Wide <p>Colors: Black or Brown</p> <p>Closure: Velcro or Lace</p> <p>Cover: Fine Leather</p>	<p>Specifications: Narrow width dress shoe. Top stitching creates shallow toe box.</p>

New Yorker



Sizes:

- 7 - 14

Widths:



- Wide
- X - Wide
- XX - Wide



Colors: Black

Closure: Lace

Cover: Napa Leather

Specifications: Same shoe pattern as Zennon, Leather, Suede, and Boot. Dress Shoe with thick toe box for added protection.

<p>Custom Diabetic Inserts</p>		
	<p>Sizes:</p> <ul style="list-style-type: none"> • 4 - 17 <p>Widths:</p> <ul style="list-style-type: none"> • Medium • Wide • X - Wide • XX - Wide <p>Options: Custom or Heat Molded</p>	<p>Custom 3-Layer Inserts Materials: Top layer is fine celled Pro-Cell Eva, the insert will naturally heat-mold to the foot without abrasion. The center layer is Pron Foam that will not bottom out and continue to provide a comfortable cushion. The bottom layer is a firmer Pro-cell Plus Eva designed to hold the shape of the insole.</p> <p>Specifications: Designed specifically for patient's foot and selected shoe. Insert works together with shoe to support the foot and arch, as well as protecting the foot's sole. Insert is specifically designed to create less friction on the foot and fits inside the shoes.</p>
<p>Gauntlet</p>		
	<p>Sizes:</p> <ul style="list-style-type: none"> • 4 - 17 <p>Widths:</p> <ul style="list-style-type: none"> • Medium • Wide • X - Wide • XX - Wide <p>Closure: Velcro</p>	<p>Foot Gauntlet</p> <p><i>For night-time relief</i></p> <ul style="list-style-type: none"> • Increases circulation • Captures body heat • Relieves pain • Provides support <p>Open heel to aid in keeping the foot dry</p>

<p>Diabetic Seamless Socks</p>		
	<p>Sizes:</p> <ul style="list-style-type: none"> • Small • Medium • Large • X - Large <p>Colors: White</p>	<p>Material: Lightweight soft cotton</p>
<p>Ritz Stick & Impression Box</p>		
	<p>Ritz Stick & Impression Box</p>	<p>Fitter uses impression box to obtain a 3-D impression of the patient's feet. The custom inserts are created from this impression. Obtaining a deep impression of the patient's feet is very important for the lab to adjust the inserts for any foot deformities. Impression box and adequate foot measurements using a Ritz stick are vital to good fitting shoes and a happy customer.</p>